


# Salgado™ HOP Driver

## INSTRUCTIONS FOR USE

Catalog # **CR2000** and **CR2005**

### Step 1

Use the malleable applicator to position the Pin over the bleeder. Push the Pin into the bone using a finger or a pair of forceps.



### Step 2


Hemorrhage Occluder Pin is inserted into the bone and applicator is removed.



### Step 3

Rotate the Salgado HOP Driver so the distal end is positioned perpendicular to the Pin head.

Gently tap downward on the top of the driver to fully seat the Pin in the sacrum.



### IMPORTANT:

- Salgado HOP Driver is supplied non-sterile.
- Salgado HOP Driver is intended for reuse. Do not use any part that appears to be defective.
- Surgin disclaims any and all responsibility and liability concerning the use of this product or any component herein, if it is reused, abused, misused, damaged or used contrary to the instructions contained herein.



10040 E. Happy Valley Rd. #415  
Scottsdale, AZ 85255 USA  
Tel: (714) 832-6300  
Fax: (714) 832-0300  
info@surgin.com  
hemoccluserpin.com

## **CLEANING INSTRUCTIONS**

### **RINSING**

Immediately after surgery, rinse instrument under warm (not hot) running water. Rinse should remove all blood, bodily fluids and tissue.

### **CLEANING**

If not done immediately after rinsing, instrument should be submerged in a solution of water and neutral PH (7) detergent.

#### **ULTRASONIC CLEANING**

- Instrument should be processed in a cleaner for the full recommended cycle time - usually 5 to 10 minutes.
- Place instrument in the ultrasonic cleaner. Instrument has to be fully submerged.
- Do not place dissimilar metals (stainless, copper, chrome plated, etc.) in the same cleaning cycle.
- Change solution frequently - at least as often as manufacturer recommends.
- Rinse instrument after ultrasonic cleaning with water to remove ultrasonic cleaning solution.

#### **AUTOMATIC WASHER STERILIZER**

Follow manufacturer's recommendations.

#### **MANUAL CLEANING**

Most instrument manufacturers recommend ultrasonic cleaning as the best and most effective way to clean surgical instruments. If ultrasonic cleaning is not available, observe the following steps:

- Use stiff plastic cleaning brushes (nylon, etc.).
- Do not use steel wood or wire brushes.
- Use only neutral PH (7) detergents because if not rinsed off properly, low PH detergents will cause breakdown of protective surface and black staining.
- High PH detergent will cause surface deposit of brown stain, which will also interfere with smooth operation of the instrument.
- Make sure all instrument surfaces are visibly clean and free from stains and tissue.
- After scrubbing, rinse instruments thoroughly under running water.

### **AFTER CLEANING**

If instrument is to be stored, let it air dry and store in a clean and dry environment.

## **STERILIZATION INSTRUCTIONS**

Surgin recommends the following AAMI guidelines for steam sterilization.

<b><u>If sterilized with porous materials:</u></b>	<b><u>If sterilized with non-porous materials:</u></b>
<b>Standard Cycle</b>	<b>Standard Cycle</b>
Sterilization time: 10 minutes	Sterilization time: 3 minutes
Temperature: 132° C (270° F)	Temperature: 132° C (270° F)
<b>Prevacuum Cycle</b>	<b>Prevacuum Cycle</b>
Sterilization time: 4 minutes	Sterilization time: 3 minutes
Temperature: 132° C (270° F)	Temperature: 132° C (270° F)